

Brucine (2,3-Dimethoxystrychnidin-10-one)

Brucine (2,3-Dimethoxystrychnidin-10-one) is an alkaloid, structurally related to strychnine, but less toxic. Figure 1 shows the ^1H NMR spectrum of a 250 mM Brucine sample in CDCl_3 measured in a single scan taking 10 seconds to acquire.

1D Proton spectrum




Brucine

Solvent = CDCl_3
Concentration = 250 mM
Frequency = 90 MHz

1D Proton

Number of scans = 1
Repetition time = 10 s
Pulse angle = 90°
Total experimental time = 10 s



10 seconds




Figure 1: ^1H NMR spectrum of a 250 mM Brucine sample in CDCl_3 measured on a Spinsolve 90 MHz system in a single scan.

1D Carbon spectrum

Figure 2 shows the ^{13}C NMR spectrum of 250 mM Brucine in CDCl_3 acquired using NOE polarization transfer from ^1H to ^{13}C and ^1H decoupling. The 1D Carbon experiment using NOE is sensitive to all ^{13}C nuclei in the sample. It clearly resolves all the expected resonances.




Brucine

Solvent = CDCl_3
Concentration = 250 mM
Frequency = 23 MHz

1D Carbon

Number of scans = 2048
Repetition time = 3.5 s
Pulse angle = 45°
Total experimental time = 120 min



120 min




Figure 2: ^{13}C NMR spectrum of a 250 mM Brucine sample in CDCl_3 measured on a Spinsolve 90 MHz system in 120 minutes.

2D COSY spectrum

The 2D COSY experiment allows one to identify coupled ^1H nuclei as they generate cross peaks out of the diagonal of the 2D data set. In Figure 2 a large number of cross peaks can be nicely observed. For example, the protons at position 6 and 11 (light green) couple with each other. Furthermore, proton 19 couples with proton 10 (light blue), 12 (orange) and 20 (pink). In addition, the couplings between protons 8 and 9 (dark blue) as well as the couplings of protons 8 and 9 with protons 14 and 15 (dark green) can be nicely observed.





Brucine

Solvent = CDCl_3
Concentration = 250 mM
Frequency = 90 MHz

COSY

Number of scans = 1
Total experimental time = 13 min



Brucine

Solvent = CDCl_3
Concentration = 250 mM
Frequency = 90 MHz

COSY

Number of scans = 1
Total experimental time = 13 min





Figure 3: ^1H 2D COSY experiment of a 250 mM Brucine sample in CDCl_3 acquired in 13 minutes on a Spinsolve 90 MHz system (top); zoom into the aliphatic region (0.5-5.0 ppm) of the ^1H 2D COSY spectrum which underlines the superb resolution.

2D JRES spectrum

This experiment is useful to identify the chemical groups generating a single line for each group by collapsing the J-coupling along the direct direction. The multiplets are generated along the vertical direction.




Figure 4: Homonuclear J-resolved (JRES) spectrum of 250 mM Brucine in CDCl₃ on a Spinsolve 90 MHz.

2D HSQC-ME

The HSQC is a powerful sequence widely used to correlate the ¹H with the one-bond coupled ¹³C nuclei. The Spinsolve is equipped with a multiplicity edited version (HSQC-ME) of this method. It provides the editing power of the DEPT-135 sequence, which is useful to differentiate between the signals of the CH₂ groups (blue) from the CH and CH₃ groups (red). Figure 5 shows the HSQC-ME spectrum of a 250 mM Brucine sample in CDCl₃ acquired in 2 minutes. The measurement time was optimized applying NUS (non uniform sampling).




Figure 5: HSQC-ME spectrum of a 250 mM Brucine sample in CDCl₃ showing the correlation between the ¹H (horizontal) and ¹³C (vertical) signals.

2D HMBC

To obtain long-range ^1H - ^{13}C correlations through two or three bond couplings, the Heteronuclear Multiple Bond Correlation (HMBC) experiment can be used. Figure 6 shows the long-range correlation of proton 8 with carbons 2, 3, 5, 7, 9 and 17 (the sequence shows the correlation with quaternary carbons, too).




Figure 6: HMBC spectrum of a 250 mM Brucine sample in CDCl_3 showing the long-range couplings between ^1H and ^{13}C nuclei.