Spinsolve 80

Brucine (2,3-Dimethoxystrychnidin-10-one)

Brucine (2,3-Dimethoxystrychnidin-10-one) is an alkaloid, structurally related to strychnine, but less toxic. Figure 1 shows the ¹H NMR spectrum of a 250 mM Brucine sample in CDCI₃ measured in a single scan taking 10 seconds to acquire.




Figure 1: 1H NMR spectrum of a 250 mM Brucine sample in CDCI₃ measured on a Spinsolve 80 MHz system in a single scan.

¹³C Spectrum

Figure 3 shows the 13 C NMR spectrum of 250 mM Brucine in CDCl₃ acquired using NOE polarization transfer from 1 H to 13 C and 1 H decoupling. The 1D Carbon experiment using NOE is sensitive to all 13 C nuclei in the sample. It clearly resolves all the expected resonances.




Figure 3: 13C NMR spectrum of a 250 mM Brucine sample in CDCl₃ measured on a Spinsolve 80 MHz system in 102 minutes.



2D COSY

Spinsolve 80

magritek

The 2D COSY experiment allows one to identify coupled ¹H nuclei as they generate cross peaks out of the diagonal of the 2D data set. In Figure 2 a large number of cross peaks can be nicely observed. For example, the protons at position 6 and 11 (light green) couple with each other. Furthermore, proton 19 couples with proton 10 (light blue), 12 (orange) and 20 (pink). In addition, the couplings between protons 8 and 9 (dark blue) as well as the couplings of protons 8 and 9 with protons 14 and 15 (dark green) can be nicely observed.




Figure 2: ¹H 2D COSY experiment of a 250 mM Brucine sample in CDCl₃ acquired in 17 minutes on a Spinsolve 80 MHz system.

2D HSQC-ME

The HSQC is a powerful sequence widely used to correlate ¹H with the one-bond coupled ¹³C nuclei. The Spinsolve is equipped with a multiplicity edited version (HSQC-ME) of this method. It provides the editing power of the DEPT-135 sequence, which is useful to differentiate the signals of CH₂ groups (blue) from CH and CH₃ groups (red). Figure 4 shows the HSQC-ME spectrum of a 250 mM Brucine in CDCl₃ acquired in 4 minutes. The measurement time was optimized applying NUS (non uniform sampling).




Figure 4: HSQC-ME spectrum of a 250 mM Brucine sample in CDCl₃ showing the correlation between the ¹H (horizontal) and ¹³C (vertical) signals.

Spinsolve 80

2D HMBC

To obtain long-range ${}^{1}H^{-13}C$ correlations through two or three bond couplings, the Heteronuclear Multiple Bond Correlation (HMBC) experiment can be used. Figure 5 shows the HMBC spectrum of a 250 mM Brucine sample measured in 34 minutes on our Spinsolve 80 MHz. As an example, the long-range correlations of proton 8 with carbons 17 (light green), 9 (dark green), 7 (dark blue), 5 (light blue), 2 (orange) and 3 (pink) are marked with circles. The experiment shows the correlation with quaternary carbons, too.




Figure 5: HMBC spectrum of a 250 mM Brucine sample in CDCl₃ showing the long-range couplings between ¹H and ¹³C nuclei.

