Spinsolve 80 ULTRA

Erythromycin A

Erythromycin A is a broad-spectrum, macrolide antibiotic on the World Health Organization's List of Essential Medicines. It is often used for the treatment of a wide variety of bacterial infections including respiratory tract infections, skin infections, pelvic inflammatory disease. Figure 1 shows the 1H NMR spectrum of a 90 mM Erythromycin A sample in CDCl, measured in a single scan taking 10 seconds to acquire.




Figure 1: 1H NMR spectrum of a 90 mM Erythromycin A sample in CDCI, measured on a Spinsolve 80 MHz ULTRA system in a single scan.

¹³C Spectrum

Figure 3 shows the ¹³C NMR spectrum of 90 mM Erythromycin A in CDCl₃ acquired using NOE polarization transfer from ¹H to ¹³C and ¹H decoupling. The 1D Carbon experiment using NOE is sensitive to all ¹³C nuclei in the sample. It clearly resolves all the expected resonances.




Figure 2: ¹³C NMR spectrum of a 90 mM Erythromycin A sample in CDCl₃ measured on a Spinsolve 80 MHz ULTRA system in 13 h 39 minutes.



The 2D COSY experiment allows one to identify coupled ¹H nuclei as they generate cross peaks out of the diagonal of the 2D data set. In Figure 3 a large number of cross peaks can be nicely observed. For example, the proton at position 1" couples to protons 2" (orange); and the proton 1' couples with proton 2' (light green). The proton 3 couples with proton 2 (dark green) and proton 4 (red). The proton 4 also couples to proton 5 (blue). In addition, the coupling between protons 14 and 15 (pink) can be observed.




Figure 3: 'H 2D COSY experiment of a 90 mM Erythromycin A sample in CDCI, acquired in 18 minutes on a Spinsolve 80 MHz ULTRA system.

2D HSQC-ME

The HSQC is a powerful sequence widely used to correlate ¹H with the one-bond coupled ¹³C nuclei. The Spinsolve is equipped with a multiplicity edited version (HSQC-ME) of this method. It provides the editing power of the DEPT-135 sequence, which is useful to differentiate the signals of CH₂ groups (blue) from CH and CH₃ groups (red). Figure 4 shows the HSQC-ME spectrum of a 90 mM Erythromycin A in CDCl₃ acquired in 137 minutes.




Figure 4: HSQC-ME spectrum of a 90 mM Erythromycin A sample in CDCl₃ showing the correlation between the ¹H (horizontal) and ¹³C



Spinsolve 80 ULTRA

2D HMBC

To obtain long-range $^1\text{H}^{-13}\text{C}$ correlations through two or three bond couplings, the Heteronuclear Multiple Bond Correlation (HMBC) experiment can be used. Figure 5 shows the HMBC spectrum of a 90 mM Erythromycin A sample in CDCl₃ measured in 274 minutes on our Spinsolve 80 MHz ULTRA. As an example, the long-range correlations of quaternary carbon 9 with protons 11 (orange), 10 (blue) and 7 (light green) are marked with circles.




Figure 5: HMBC spectrum of a 90 mM Erythromycin A sample in CDCl₃ showing the long-range couplings between ¹H and ¹³C nuclei.

