

## Erythromycin A

Erythromycin A is a broad-spectrum, macrolide antibiotic on the World Health Organization's List of Essential Medicines. It is often used for the treatment of a wide variety of bacterial infections including respiratory tract infections, skin infections, pelvic inflammatory disease. Figure 1 shows the <sup>1</sup>H NMR spectrum of a 90 mM Erythromycin A sample in CDCl<sub>3</sub> measured in a single scan taking 10 seconds to acquire.




Figure 1: <sup>1</sup>H NMR spectrum of a 90 mM Erythromycin A sample in CDCl<sub>3</sub> measured on a Spinsolve 80 MHz ULTRA system in a single scan.

## <sup>13</sup>C Spectrum

Figure 3 shows the <sup>13</sup>C NMR spectrum of 90 mM Erythromycin A in CDCl<sub>3</sub> acquired using NOE polarization transfer from <sup>1</sup>H to <sup>13</sup>C and <sup>1</sup>H decoupling. The 1D Carbon experiment using NOE is sensitive to all <sup>13</sup>C nuclei in the sample. It clearly resolves all the expected resonances.




Figure 2: <sup>13</sup>C NMR spectrum of a 90 mM Erythromycin A sample in CDCl<sub>3</sub> measured on a Spinsolve 80 MHz ULTRA system in 13 h 39 minutes.

## 2D COSY

## Spinsolve 80<sup>ULTRA</sup>

The 2D COSY experiment allows one to identify coupled  $^1\text{H}$  nuclei as they generate cross peaks out of the diagonal of the 2D data set. In Figure 3 a large number of cross peaks can be nicely observed. For example, the proton at position 1'' couples to protons 2'' (orange); and the proton 1' couples with proton 2' (light green). The proton 3 couples with proton 2 (dark green) and proton 4 (red). The proton 4 also couples to proton 5 (blue). In addition, the coupling between protons 14 and 15 (pink) can be observed.



### Erythromycin A

Solvent =  $\text{CDCl}_3$   
Concentration = 90 mM  
Frequency = 80 MHz

### COSY

Number of scans = 1  
Number of steps = 512  
Total experimental time = 18 min





Figure 3:  $^1\text{H}$  2D COSY experiment of a 90 mM Erythromycin A sample in  $\text{CDCl}_3$  acquired in 18 minutes on a Spinsolve 80 MHz ULTRA system.

## 2D J-Res

The 2D J-Res experiment is useful to identify the chemical groups generating a single line for each group by collapsing the J-coupling along the direct direction. The multiplets are generated along the vertical direction.



### Erythromycin A

Solvent =  $\text{CDCl}_3$   
Concentration = 90 mM  
Frequency  $^1\text{H}$  = 80 MHz

### JRES

Number of scans = 4  
Total experimental time = 43 min





Figure 4:  $^1\text{H}$  2D J-Res experiment of a 90 mM Erythromycin A sample in  $\text{CDCl}_3$  acquired in 43 minutes on a Spinsolve 80 MHz ULTRA system.



## 2D HSQC-ME

## Spinsolve 80 <sup>ULTRA</sup>

The HSQC is a powerful sequence widely used to correlate  $^1\text{H}$  with the one-bond coupled  $^{13}\text{C}$  nuclei. The Spinsolve is equipped with a multiplicity edited version (HSQC-ME) of this method. It provides the editing power of the DEPT-135 sequence, which is useful to differentiate the signals of  $\text{CH}_2$  groups (blue) from  $\text{CH}$  and  $\text{CH}_3$  groups (red). Figure 5 shows the HSQC-ME spectrum of a 90 mM Erythromycin A in  $\text{CDCl}_3$  acquired in 137 minutes.




### Erythromycin A

Solvent =  $\text{CDCl}_3$   
Concentration = 90 mM  
Frequency  $^1\text{H}$  = 80 MHz

### HSQC-ME

Number of scans = 16  
Repetition time = 1 s  
Number of steps = 512  
Total experimental time = 137 min



Red =  $\text{CH}$  and  $\text{CH}_3$

Blue =  $\text{CH}_2$

\* residual solvent




Figure 5: HSQC-ME spectrum of a 90 mM Erythromycin A sample in  $\text{CDCl}_3$  showing the correlation between  $^1\text{H}$  and  $^{13}\text{C}$  signals.

## 2D HMBC

To obtain long-range  $^1\text{H}$  -  $^{13}\text{C}$  correlations through two or three bond couplings, the Heteronuclear Multiple Bond Correlation (HMBC) experiment can be used. Figure 6 shows the HMBC spectrum of a 90 mM Erythromycin A sample in  $\text{CDCl}_3$  measured in 274 minutes on our Spinsolve 80 MHz ULTRA. As an example, the long-range correlations of quaternary carbon 9 with protons 11 (orange), 10 (blue) and 7 (light green) are marked with circles.




### Erythromycin A

Solvent =  $\text{CDCl}_3$   
Concentration = 90 mM  
Frequency  $^1\text{H}$  = 80 MHz

### HMBC

Number of scans = 64  
Repetition time = 1 s  
Number of steps = 256  
Total experimental time = 274 min



\* residual solvent




Figure 6: HMBC spectrum of a 90 mM Erythromycin A sample in  $\text{CDCl}_3$  showing the long-range couplings between  $^1\text{H}$  and  $^{13}\text{C}$  nuclei.

