

## Brucine (2,3-Dimethoxystrychnidin-10-one)

Brucine (2,3-Dimethoxystrychnidin-10-one) is an alkaloid, structurally related to strychnine, but less toxic. Figure 1 shows the <sup>1</sup>H NMR spectrum of a 50 mM brucine sample in CDCl<sub>3</sub> measured in a single scan taking 15 seconds to acquire.




Figure 1: <sup>1</sup>H NMR spectrum of a 50 mM brucine sample in CDCl<sub>3</sub> measured on a Spinsolve 80 MHz system in a single scan.

## 2D COSY

The 2D COSY experiment allows one to identify coupled <sup>1</sup>H nuclei as they generate cross peaks out of the diagonal of the 2D data set. In Figure 2 a large number of cross peaks are well identifiable. For example, the protons at position 6 and 11 (light green) couple with each other. Furthermore, proton 19 couples with proton 10 (light blue), 12 (orange) and 20 (pink). In addition, the couplings between protons 8 and 9 (dark blue) as well as the couplings of protons 8 and 9 with protons 14 and 15 (dark green) can be observed nicely.




Figure 2: <sup>1</sup>H 2D COSY experiment of a 50 mM brucine sample in CDCl<sub>3</sub> acquired in 17 minutes on a Spinsolve 80 MHz system.

### <sup>13</sup>C Spectrum

Figure 3 shows the <sup>13</sup>C NMR spectrum of 1 M brucine in CDCl<sub>3</sub> acquired using NOE polarization transfer from <sup>1</sup>H to <sup>13</sup>C and <sup>1</sup>H decoupling. The 1D carbon experiment using NOE is sensitive to all <sup>13</sup>C nuclei in the sample. It clearly resolves all the expected resonances.




Figure 3: <sup>13</sup>C NMR spectrum of a 1 M brucine sample in CDCl<sub>3</sub> measured on a Spinsolve 80 MHz system in 6.5 minutes.

### 2D HSQC-ME and 2D HMBC

The HSQC and HMBC are powerful sequences widely used to correlate the one-bond coupled <sup>1</sup>H-<sup>13</sup>C nuclei and long range <sup>1</sup>H-<sup>13</sup>C, respectively. The Spinsolve is equipped with a multiplicity edited version (HSQC-ME) of this method. It provides the editing power of the DEPT-135 sequence, which is useful to differentiate between the signals of the CH<sub>2</sub> groups (blue) from the CH and CH<sub>3</sub> groups (red). Thanks to the higher sensitivity of these sequences, compared to direct detected <sup>13</sup>C experiments, samples with lower concentrations can be measured in relatively short times. Figure 4 shows the HSQC-ME (left) and HMBC (right) spectra of a 250 mM brucine sample in CDCl<sub>3</sub> acquired in 4 minutes and 34 minutes, respectively. The HSQC spectrum resolves all coupled <sup>13</sup>C and allows one to resolve the <sup>1</sup>H signals that overlap in the 1D proton spectrum. On the other hand, the HMBC shows the long-range correlations that are useful to identify and assign the different peaks in the spectrum. As an example, we marked on the HMBC spectrum the cross peaks of proton 8 with carbons 2, 3, 5, 7, 9 and 17 (the sequence shows the correlation with quaternary carbons, too).




Figure 4: HSQC-ME (left) and HMBC (right) spectra of a 250 mM brucine sample in CDCl<sub>3</sub> showing the correlation between the <sup>1</sup>H (horizontal) and <sup>13</sup>C (vertical) signals.