

## Levofloxacin

Levofloxacin is a commonly known antibiotic drug active against a range of bacterial infections. Structurally Levofloxacin belongs to the group of fluoroquinolones. It is a chiral compound in the (-)-(S) enantiomeric form. Figure 1 shows the <sup>1</sup>H NMR spectrum of a 250 mM Levofloxacin sample in DMSO-d<sub>6</sub> measured in a single scan taking 10 seconds to acquire.

### 1D Proton spectrum




#### Levofloxacin

Solvent = DMSO-d<sub>6</sub>  
Concentration = 250 mM  
Frequency = 90 MHz

#### 1D Proton

Number of scans = 1  
Repetition time = 10 s  
Pulse angle = 90°  
Total experimental time = 10 s



10 seconds




Figure 1: <sup>1</sup>H NMR spectrum of a 250 mM Levofloxacin sample in DMSO-d<sub>6</sub> measured on a Spinsolve 90 MHz system in a single scan.

### 1D Carbon spectrum

Figure 2 shows the <sup>13</sup>C NMR spectrum of 250 mM Levofloxacin in DMSO-d<sub>6</sub> acquired using NOE polarization transfer from <sup>1</sup>H to <sup>13</sup>C and <sup>1</sup>H decoupling. The 1D Carbon experiment using NOE is sensitive to all <sup>13</sup>C nuclei in the sample. It clearly resolves all the expected resonances. The spectrum was acquired with both <sup>1</sup>H and <sup>19</sup>F decoupling.




#### Levofloxacin

Solvent = DMSO-d<sub>6</sub>  
Concentration = 250 mM  
Frequency = 90 MHz

#### 1D Carbon (19F)

Number of scans = 2048  
Repetition time = 3 s  
Pulse angle = 90°  
Total experimental time = 102 min



102 min




Figure 2: <sup>13</sup>C NMR spectrum of a 250 mM Levofloxacin sample in DMSO-d<sub>6</sub> measured on a Spinsolve 90 MHz system in 102 minutes.

## 2D COSY spectrum

The 2D COSY experiment allows one to identify coupled  $^1\text{H}$  nuclei as they generate cross peaks out of the diagonal of the 2D data set. In Figure 3 a large number of cross peaks can be observed nicely. For example, the protons at position 26 and 19 (light blue) couple with each other. Furthermore, protons 12 and 13 couple with protons 14 and 15 (orange). In addition, the couplings between proton 27 and 5 (dark blue), 19 (dark green) and 26 (light green) can be observed nicely.




Figure 3:  $^1\text{H}$  2D COSY experiment of a 250 mM Levofloxacin sample in DMSO- $d_6$  acquired in 8.5 minutes on a Spinsolve 90 MHz system.

## 2D $^1\text{H}$ - $^{13}\text{C}$ HSQC-ME

The HSQC is a powerful sequence widely used to correlate the  $^1\text{H}$  with the one-bond coupled  $^{13}\text{C}$  nuclei. The Spinsolve is equipped with a multiplicity edited version (HSQC-ME) of this method. It provides the editing power of the DEPT-135 sequence, which is useful to differentiate between the signals of the  $\text{CH}_2$  groups (blue) from the  $\text{CH}$  and  $\text{CH}_3$  groups (red). Figure 4 shows the HSQC-ME spectrum of a 250 mM Levofloxacin sample in DMSO- $d_6$  acquired in 4 minutes. The measurement time was optimized applying NUS (non uniform sampling).






Figure 4: HSQC-ME spectrum of a 250 mM Levofloxacin sample in DMSO- $d_6$  showing the correlation between the  $^1\text{H}$  (horizontal) and  $^{13}\text{C}$  (vertical) signals.

## 2D $^1\text{H}$ - $^{13}\text{C}$ HMBC

To obtain long-range  $^1\text{H}$ - $^{13}\text{C}$  correlations through two or three bond couplings, the Heteronuclear Multiple Bond Correlation (HMBC) experiment can be used. Figure 5 shows the long-range correlation of proton 5 with carbons 1, 4, 6, 18 and 23 (the sequence shows the correlation with quaternary carbons, too).



## 1D Fluorine spectrum



## 2D $^1\text{H}$ - $^{15}\text{N}$ HMBC

To obtain long-range  $^1\text{H}$ - $^{15}\text{N}$  correlations through two or three bond couplings, the Heteronuclear Multiple Bond Correlation (HMBC) experiment can be used. Figure 7 shows the full spectral range and depicts nicely all long-range correlations of the three nitrogen nuclei 3, 11 and 16. Here, nitrogen 3 couples with the protons at positions 5 and 26, nitrogen 11 couples with protons 12 and 13, as well as nitrogen 16 couples with protons 14,15 and 21, respectively. To note is the great SNR on the  $^{15}\text{N}$  trace that indicates even reduced total measurement times to the 2 h 17 min shown.




Figure 7:  $^{15}\text{N}$ -HMBC spectrum of a 250 mM Levofloxacin sample in DMSO- $d_6$  showing the long-range couplings between  $^1\text{H}$  and  $^{15}\text{N}$  nuclei.