

Bringing NMR to the Battery Lab

Practical Electrolyte Analysis with Benchtop NMR Systems



Lithium-ion batteries (LIBs) have become indispensable in modern portable electronics and are now central to the rapid expansion of electric mobility and advanced energy storage systems. Their performance, safety, and longevity depend strongly on the chemical composition and stability of the liquid electrolyte, which typically consists of a lithium salt such as LiPF_6 dissolved in mixtures of organic carbonates. Serving as the charge-transport medium between electrodes, the electrolyte governs key battery characteristics including ion mobility, rate capability, and overall energy release [1]. Understanding and monitoring electrolyte properties is therefore essential in fields ranging from fundamental materials research to industrial production and quality control.

Nuclear magnetic resonance (NMR) spectroscopy has emerged as one of the most informative analytical techniques for studying LIB electrolytes. It provides direct molecular-level insights into ion coordination, speciation, solvent interactions, degradation pathways, and transport phenomena. Numerous studies have demonstrated the value of NMR for quantifying electrolyte composition, tracking aging processes, and probing lithium-ion mobility under different conditions.

Modern Spinsolve benchtop NMR spectrometers make advanced NMR capabilities accessible in the laboratory environment with an unmatched simplicity & cost. It offers high-quality relevant data in a compact, cryogen-free platform. The system is engineered for simplicity and robustness, enabling reliable measurements directly at the point of use. Neat electrolyte samples can be analyzed without dilution or deuterated solvents, thanks to the integrated hardware lock. Based on a go-without-saying robustness, Spinsolve also provides fully quantitative results with the use of a unique (per nucleus) internal or external

standard saving cumbersome calibration time & standard cost. Spinsolve supports automatic nucleus switching among key nuclei: ^1H , ^7Li , ^{11}B , ^{13}C , ^{19}F , ^{23}Na , and ^{31}P , eliminating manual tuning and matching steps.

With a single unit, Spinsolve combines identification, quantitative & mobility results on solvents, additives, conducting salts and their degradation products. Spinsolve consolidates workflows that traditionally require multiple analytical techniques. Its speed, operational ease, and minimal maintenance demands make it an attractive platform for routine battery research, production control, and diagnostic applications.

In this Application Note, we highlight how Spinsolve can be used to rapidly and effectively characterize lithium-ion battery electrolytes, demonstrating its value as a versatile tool for both research laboratories and industrial environments.

• Solvent composition analysis

LIB electrolytes usually consist of ethylene carbonate (EC) and other carbonate solvents. EC with high dielectric constant facilitates disassociation of (most often LiPF_6) salts to give free ions. It is also involved in the formation of the solid electrolyte interface (SEI) to protect the anode. But the viscous nature of EC requires blending with other carbonate solvents, commonly diethylene carbonate (DEC), dimethyl




Figure 1. Quantification of EC, DMC, and DEC in a ternary electrolyte mixture by ^1H NMR. The comparison against gravimetric quantification shows excellent agreement (table top right). The top left plot shows an overlay of 60 spectra measured to assess the repeatability of the NMR analysis. The relative standard deviation of the 60 measurements is less than 0.1% for all three components.

carbonate (DMC), ethylmethyl carbonate (EMC), or propylene carbonate (PC) to increase the ion mobility and enhance temperature performance. The composition of solvents needs to be carefully tailored for a good balance and accurately controlled to ensure optimum performance. Figure 1 shows a single-scan ^1H spectrum of a mixture of DMC, DEC and EC. The signals of all three solvents can be baseline-resolved. The characteristic chemical shifts and splitting pattern allows one to analyze solvent composition of any unknown electrolytes. Based on the integrals, the ratio of each solvent can be easily calculated with high precision ($\text{RSD} < 0.1\%$) and is in excellent agreement with gravimetric results. The spectrum was recorded in 7 seconds without any sample preparation prior to NMR measurement, offering a convenient way for solvent composition analysis compared to high performance liquid chromatography (HPLC) or gas chromatography (GC).

- **Fast detection and quantification of additives**

The high sensitivity of ^1H NMR also allows detecting and quantifying lower concentration of additives, such as vinylene carbonate (VC) as a common additive in electrolytes to prolong battery lifetime. To highlight the linearity and repeatability of the Spinsolve system, a series of DMC/DEC/EC electrolytes with 0-3 % of VC were measured. VC contributes to a characteristic singlet at 7.2 ppm. The Spinsolve spectrometer has a fixed response factor for individual nuclei, the signal integral is independent of parameter settings and is directly proportional to molar concentration—not only for components within the same mixture, but also for samples measured in separate tubes. Figure 2 shows a series of electrolytes spiked with various amounts of VC. The integral of the VC signal shows an excellent linear correlation with the VC content. Because of the high sensitivity, a one-minute ^1H NMR measurement delivered quantification results with




Figure 2. quantification of VC spiked in electrolytes by ^1H NMR.

high repeatability for an electrolyte containing 1 % VC (1.000 ± 0.0025 %, based on 10 consecutive measurements). A recent study [2] reported quantification of 1-5 % VC in electrolytes by Spinsolve NMR spectrometer to investigate the amount of VC involved in SEI formation. Other widely used additives, like ethylene sulfate (DTD), 1,3,6-hexanetricarbonitrile (HTN), and fluoroethylene carbonate (FEC) can likewise be identified and quantified by using ^1H NMR spectroscopy, as shown in ^1H spectra of some other complex electrolytes in Figure 3.




Figure 3. ^1H spectra of artificial electrolyte (in red) containing DTD and HTN in EMC and electrolyte (in blue) containing VC and FEC in EC and EMC.

- **Quantification of conducting salts**

Quantitative analysis using techniques such as HPLC, GC, or high-field NMR typically requires frequent calibration with reference standards. On a Spinsolve spectrometer, the response factor for each nucleus is determined by the factory hardware settings and remains highly stable over time due to the exceptional stability of the magnet and electronics. As a result, repeated calibration is generally unnecessary. In addition, the response factor is minimally affected by sample properties, as evidenced by the excellent linearity of Li^+ over a wide dynamic range shown in Figure 4.

As the initial step, a single-point calibration using an external standard is sufficient to determine the response factor, which can be used to calculate the concentration of unknown samples based on the integral. In this example, the Li^+ response factor was first established with a solution of lithium chloride (0.5 M in water). The Li^+ concentration in the same electrolytes with and without aging treatment were easily derived from the integrals. The lower Li^+ concentration in aged electrolyte can be attributed to lithium salt precipitation.




Figure 4. Regression curve of LiCl at different concentrations and Li^+ quantification of electrolytes using a single-point calibration.

This quantification approach can be used not only to quantify the Li^+ cation concentration but also the anion concentration. This is particularly important for multi-salt electrolytes designed to enhance the electrochemical and thermal stability of advanced batteries [3]. Because the anions in these systems typically contain fluorine, ^{19}F NMR spectroscopy is especially powerful due to its high sensitivity and broad chemical-shift dispersion.




Figure 5. ^{19}F spectrum of a multi-salt electrolyte showing signals from conducting salts (marked regions) and their related hydrolysis products (minor signals outside of marked regions).

This wide dispersion allows ^{19}F NMR to resolve multiple fluorine-containing species within the same spectrum, making it well suited for characterizing complex multi-salt electrolyte formulations. Note that the signals are about 200 ppm apart, uniform excitation of all signals on high-field NMR can be difficult, even if the transmitter frequency is carefully chosen at the center. The Spinsolve Multi X switching technology combined with the efficient probe design enables a uniform excitation profile over a large chemical shift range enabling accurate quantification of signals that are even more than 200 ppm apart from each other without the need of correction or additional measurements. Figure 5 shows the ^{19}F spectrum of an artificial electrolyte where all three lithium-conducting salts are clearly resolved in a distinct chemical-shift region, along with their corresponding hydrolysis products. As only a single response factor is required for quantification, the integral values of the different anions can be directly converted into molar concentrations that can be further converted in mass concentration once the structures are identified.

- **Hydrofluoric acid (HF) quantification**

Quantifying HF in electrolytes is essential because a small amount of HF can form through oxidation processes and severely affect the chemical stability, safety, and performance of electrochemical systems. HF is highly reactive and corrosive, and it can degrade electrode materials, damage protective surface layers, and trigger unwanted side reactions. This leads to accelerated aging, reduced efficiency, and potential safety risks. By monitoring HF levels as part of routine quality control, manufacturers ensure material integrity, consistent product performance, and long-term reliability.




Figure 6. HF quantification in fresh and aged electrolytes by ^{19}F measurement. The signal of HF appears as a doublet in a fresh electrolyte because of ^1H - ^{19}F coupling but collapses to a broad singlet due to exchange at higher water content.

Figure 6 shows ^{19}F spectra of a commercial electrolyte, a doublet at -74.3 ppm from PF_6 anion dominates the spectrum. After more scans, however, the signal of HF gets also resolved at 190.1 ppm. In a fresh electrolyte, the HF signal appears as a doublet because HF remains intact at very low water content, and the coupling to the bonded proton causes splitting of the fluorine signal. In the presence of larger water content, the doublet collapses into a singlet because HF undergoes rapid proton exchange with water and related species, averaging out the ^{19}F - ^1H coupling. Based on its integral, the HF content in the fresh

electrolyte was calculated to be 40.2 ppm. The signal-to-noise ratio (SNR) of the HF signal is 12, defining a limit of detection (LOD – concentration for SNR=3), of 10 ppm for a 8 minutes measurement. Considering that the SNR of the NMR spectrum is proportional to the square root of measurement time, the LOD can be further reduced by increasing the measurement time.

- Investigation of electrolyte degradation

Detecting and quantifying HF and other hydrolysis products formed from the reaction between the conducting salt and moisture in air or residual water in the electrolyte is important for understanding electrolyte aging. Figure 7 shows a ^{19}F spectrum of an aged electrolyte acquired in a single-scan (5 seconds). In the spectrum we can clearly identify the PF_6^- anion and its related hydrolysis products. Signals were assigned following previous study [4,5].




Figure 7. ^{19}F spectrum of an aged electrolyte acquired in just 5 seconds. In this short measurement time different LiPF_6 hydrolysis products were detected.

An accelerated PF_6 degradation study was conducted by adding two microliters of water into 0.5 mL of a standard LiPF_6 electrolyte. Figure 8 shows the concentration as a function of time of the different hydrolysis products. The monitoring experiment acquired ^{19}F spectra during 20 hours. Based on the proposed degradation pathway [5], if water is abundant, PF_6 is hydrolyzed rapidly to OPF_2O^- and strongly increases HF concentration. OPF_2O^- further reacts to water to generate a more stable OPFO_2^{2-} and releases more HF. HF attacks the borosilicate glass tube used to contain the sample to yield BF_4^- and water. The reaction between PF_6 -derived related species (such as OPF_2O^-) and solvent molecules like DMC is chemically feasible, the resulting products typically only reach detectable levels after an extended aging period due to the extremely slow reaction rate [5].




Figure 8. Accelerated hydrolysis of electrolytes monitored by ^{19}F NMR.

- **Characterization of PVDF**

Polyvinylidene fluoride (PVDF) based electrolytes are of great potential for the development of solid-state batteries [7]. It functions also as separator component in LIB [8] and as binder for electrode materials [9]. The ^{19}F spectrum of PVDF provides valuable structural information. The highest signal between -90 and -95 ppm corresponds to the head-to-tail connection, while the head-to-head connection form is also visible between -100 and -120 ppm [10]. The ratio between the two forms can be simply calculated as the ratio of the respective signal integrals. As the different ^{19}F signals are excited uniformly, no complex correction is required, as explained above. The fact that the ratio of the signals at around -72 ppm (CF_3 group) and -180 ppm (CF group) is exactly 3:1, suggests the presence of branching structures from the hexafluoropropylene (HFP) monomer. The signals of the HFP- CF_2 group, which should exhibit an integral of 2, overlap with the VF_2 - CF_2 signals from -90 to -120 ppm. By subtracting the HFP- CF_2 integral, the integral of VF_2 - CF_2 can be calculated as 7.5. Then the monomer ratio between VF_2 and HFP is 3.755:1. Combining the monomer ratio and the fluorine content of each monomer (VF_2 : 59.3%; HFP: 76.0%), the total fluorine content of the PVDF polymer can be also calculated as 65.7%.




Figure 9. ^{19}F NMR spectrum of PVDF-HFP copolymer in DMF. Spectral assignments: -65 to -75 ppm (HFP- CF_3 groups); -90 to -95 ppm (VDF- CF_2 head-to-tail connection); -100 to -120 ppm (HFP- CF_2 groups overlapped with VDF- CF_2 head-to-head connection); and -175 to -185 ppm (HFP-CF groups).

- Ion mobility and other performance indicators of electrolytes

Spinolve NMR spectrometers can be equipped with pulsed field gradient (PFG) capabilities to enable the measurement of molecular self-diffusion coefficients. The diffusion coefficients derived from PFG NMR can be related to ion mobility, ion-ion interaction, solvation effect, and ultimately to battery performance. Figure 10 summarizes the diffusion measurements of all three nuclei in a LiPF₆ electrolyte (1M LiPF₆ in EC/DMC 1:1, v/v). As expected from ¹⁹F and ³¹P measurements, a similar diffusion value is obtained because they diffuse together. Surprisingly, lithium ions diffuse slower than the anions although the size is smaller. This is due to the large charge density of Li⁺, which attracts solvent molecules to form a big complex [11].




Figure 10. Diffusion NMR measurement of cation and anion in electrolyte.

In addition, key performance indicators of the ion transport property such as transference number and ion conductivity can be derived from diffusion values based on the Nernst–Einstein equation [12]. Table 1 summarizes the results for three groups of electrolytes consisting of 1 M LiPF₆ in different solvent systems. The transference numbers, t_+ and t_- , show little dependence on solvent composition and ratio. EC/DMC electrolytes show higher conductivity than EC/DEC, while the conductivity of the EC/DEC/DMC electrolyte falls between the two. Between the two EC/DMC electrolytes, a higher proportion of DMC apparently further increases the ion conductivity.

It has been noted that Nernst–Einstein equation overestimates the ion conductivity [12,13]. As shown in Table 1, the conductivity calculated from NMR diffusion measurements (σ_{NMR}) is higher than the value reported by the supplier (σ_{real} , <https://solvionic.com/en/6-electrolytes>). This discrepancy appears because the Nernst–Einstein equation assumes an ideal model where lithium salts are fully dissociated and all ions diffuse independently, which, in reality, is rarely the case. For example, ion interactions lead to ion

pairs and ion clusters, especially at such high salt concentrations, effectively reducing the number of charge carriers. Nevertheless, the $\sigma_{\text{NMR}}/\sigma_{\text{real}}$ ratio remains approximately constant at about 2 for the three electrolytes, indicating that the conductivity results obtained by NMR could be a reliable comparative performance indicator. The $\sigma_{\text{NMR}} / \sigma_{\text{real}}$ ratio is linked often to the degree of ion association [14]. A lower $\sigma_{\text{NMR}} / \sigma_{\text{real}}$ ratio reflects a smaller extent of ion pairing or aggregation and weaker ion-ion or ion-solvent interactions, corresponding to a larger fraction of free ions actively contributing to charge transport. Therefore, this ratio can serve as quantitative parameter for evaluating dissociation and guiding electrolyte design.

Table 1. Summary of diffusion values and calculated performance indicators.

Solvent	Solvent ratio	D value (10 ⁻¹⁰ m ² /s)			Performance indicators			
		⁷ Li	¹⁹ F	³¹ P	t+	t-	Conductivity σ (mS/cm)	
							NMR	Real
EC/DEC	1:1	1.67	2.77	2.72	0.38	0.62	16.6	8.07
	3:7	2.00	3.00	2.99	0.40	0.60	18.7	-
EC/DMC	1:1	2.40	3.80	3.76	0.39	0.61	23.2	11.38
	3:7	3.06	4.39	4.43	0.41	0.59	27.8	-
EC/DMC/DEC	1:1:1	2.07	3.11	3.11	0.40	0.60	19.4	9.65

Conclusions

This Application Note demonstrates that Spinsolve benchtop NMR spectrometers provide a uniquely powerful and practical platform for comprehensive characterization of lithium-ion battery electrolytes. By combining high spectral quality, true quantitative measurements, fast measurement times, and exceptional operational simplicity in a compact, cryogen-free instrument, Spinsolve spectrometers enable analytical capabilities that were previously restricted to high-field NMR laboratories. The ability to measure neat electrolytes without deuterated solvents, automatically switch between multiple nuclei, and obtain accurate quantitative results without demanding calibration, greatly enhances the analysis of electrolytes—from solvent composition and additive quantification to salt concentration, degradation monitoring, and ion mobility studies.

Across all examples shown in this application note, ¹H, ⁷Li, ¹⁹F, and ³¹P spectroscopy, as well as diffusion NMR, Spinsolve spectrometers deliver the resolution, sensitivity, and long-term stability needed for both advanced research and routine quality control. Complex electrolyte formulations, multi-salt systems, and degradation pathways can be examined quickly and with high confidence, while gradient-based diffusion measurements provide deeper insight into ion transport and electrolyte performance metrics. Bringing these capabilities together in a single, easy-to-operate instrument reduces cost, minimizes sample preparation needs, and lowers operational complexity compared to conventional approaches.

As the pace of battery innovation continues to accelerate, the demand for fast, accessible, and nondestructive analytical techniques grows accordingly. Spinsolve benchtop NMR brings high-value molecular insights directly to the laboratory or production environment, enabling researchers, engineers, and manufacturers to shorten development cycles, strengthen quality assurance, and advance their understanding of electrolyte behavior. Its robustness, ease of use, and wide analytical scope make Spinsolve an indispensable tool for the next generation of battery science and industry.

References

- [1] K. Xu, "Nonaqueous Liquid Electrolytes for Lithium-Based Rechargeable Batteries," *ChemInform*, vol. 35, no. 50, Dec. 2004, doi: 10.1002/chin.200450271.
- [2] S. Azam, W. Black, H. MacLennan, A. Eldesoky, and J. R. Dahn, "Additive Screening of LFP/graphite Pouch Cells for High Temperature Cycling at 70 °C," *J. Electrochem. Soc.*, vol. 172, no. 2, p. 020536, 2025, doi: 10.1149/1945-7111/adb64d.
- [3] J. Hu *et al.*, "Rational design of multi-salt propylene carbonate-based electrolyte for enhanced long-lifespan and safety in pouch cells," *Chemical Engineering Journal*, vol. 518, p. 164639, Aug. 2025, doi: 10.1016/j.cej.2025.164639.
- [4] P. Barnes *et al.*, "A non-aqueous sodium hexafluorophosphate-based electrolyte degradation study: Formation and mitigation of hydrofluoric acid," *J. Power Sources*, vol. 447, p. 227363, 2020, doi: <https://doi.org/10.1016/j.jpowsour.2019.227363>.
- [5] S. Wiemers-Meyer, M. Winter, and S. Nowak, "Mechanistic insights into lithium ion battery electrolyte degradation-a quantitative NMR study," 2016, *Royal Society of Chemistry*. doi: 10.1039/c6cp05276b.
- [6] J. Hu *et al.*, "Rational design of multi-salt propylene carbonate-based electrolyte for enhanced long-lifespan and safety in pouch cells," *Chemical Engineering Journal*, vol. 518, p. 164639, 2025, doi: <https://doi.org/10.1016/j.cej.2025.164639>.
- [7] Y. Wu, Y. Li, Y. Wang, Q. Liu, Q. Chen, and M. Chen, "Advances and prospects of PVDF based polymer electrolytes," *Journal of Energy Chemistry*, vol. 64, pp. 62–84, 2022, doi: <https://doi.org/10.1016/j.jechem.2021.04.007>.
- [8] K. Bicy, A. B. Gueye, D. Rouxel, N. Kalarikkal, and S. Thomas, "Lithium-ion battery separators based on electrospun PVDF: A review," *Surfaces and Interfaces*, vol. 31, p. 101977, 2022, doi: <https://doi.org/10.1016/j.surfin.2022.101977>.
- [9] X. Zhong *et al.*, "Binding mechanisms of PVDF in lithium ion batteries," *Appl. Surf. Sci.*, vol. 553, p. 149564, 2021, doi: <https://doi.org/10.1016/j.apsusc.2021.149564>.
- [10] S. Banerjee, Y. Patil, O. Gimello, and B. Ameduri, "Well-defined multiblock poly (vinylidene fluoride) and block copolymers thereof: a missing piece of the architecture puzzle," *Chemical Communications*, vol. 53, no. 79, pp. 10910–10913, 2017.
- [11] Z. Cui *et al.*, "Molecular anchoring of free solvents for high-voltage and high-safety lithium metal batteries," *Nat. Commun.*, vol. 15, no. 1, p. 2033, 2024.
- [12] A. France-Lanord and J. C. Grossman, "Correlations from Ion Pairing and the Nernst-Einstein Equation," *Phys. Rev. Lett.*, vol. 122, no. 13, Apr. 2019, doi: 10.1103/PhysRevLett.122.136001.
- [13] M. Chintapalli *et al.*, "Relationship between Conductivity, Ion Diffusion, and Transference Number in Perfluoropolyether Electrolytes," *Macromolecules*, vol. 49, no. 9, pp. 3508–3515, May 2016, doi: 10.1021/acs.macromol.6b00412.

Table of common electrolyte components quantified with Spinsolve

Category	Compound	Full name	Nuclei	Comment
Carbonate solvents	EC	Ethylene carbonate	^1H , ^{13}C	-
	DMC	Dimethyl carbonate	^1H , ^{13}C	-
	DEC	Diethylene carbonate	^1H , ^{13}C	-
	EMC	Ethylmethyl carbonate	^1H , ^{13}C	-
	PC	Propylene carbonate	^1H , ^{13}C	-
Additives	VC	Vinylene carbonate	^1H	LOD* 18 ppm
	FEC	Fluoroethylene carbonate	^1H , ^{19}F	^{19}F LOD 38 ppm
	DTD	Ethylene sulfite	^1H	LOD 12 ppm
	PES	1,3-Propene sultone	^1H	
	PDO	3-Phenyl-1,4,2-dioxazol-5-one	^1H	
	HTN	1,3,6-Hexanetricarbonitrile	^1H	LOD 85 ppm
	LiDFOB	Lithium difluoro(oxalato)borate	^1H , ^7Li , ^{11}B , ^{19}F	^{19}F LOD 99 ppm
	LiDFP	Lithium difluorophosphate	^1H , ^7Li , ^{19}F , ^{31}P	^{19}F LOD 64 ppm
Conducting salts	LiPF ₆	Lithium hexafluorophosphate	^7Li , ^{19}F , ^{31}P	-
	LiFSI	Lithium bis(fluorosulfonyl)imide	^7Li , ^{19}F , ^{31}P	-
	LITFSI	Lithium bis(trifluoromethanesulfonyl)imide	^7Li , ^{13}C , ^{19}F , ^{31}P	-
	LiBF ₄	Lithium tetrafluoroborate	^7Li , ^{11}B , ^{19}F	-
Hydrolysis products	HF	Hydrofluoric acid	^{19}F	^{19}F LOD 10 ppm
	LiDFP	Lithium difluorophosphate	^1H , ^7Li , ^{19}F , ^{31}P	^{19}F LOD 64 ppm
	LiPO ₃ F	Lithium fluorophosphate	^1H , ^7Li , ^{19}F , ^{31}P	^{19}F LOD 120 ppm
Battery component	H ₂ O	Residual water	^1H	LOD 5 ppm
	PVDF-HFP	Poly(vinylidene fluoride-co-hexafluoropropylene)	^{19}F	-

*LOD values are given for 2 min measurement time. Longer measurement times result in the reduction of the LODs.

This table does not cover all components detectable with Spinsolve spectrometers. If the species you are interested in are not included, please contact sales@magritek.com for assistance.

Contact us now for a quote, to request a demo or to measure your samples

Email: sales@magritek.com

Website: www.magritek.com/contact-us

GERMANY +49 241 9278 7270

UNITED STATES +1 855 667 6835

UNITED KINGDOM +44 7468 529 615

NEW ZEALAND +64 4 477 7096

For a complete listing of our global offices and distributors visit: www.magritek.com/about-us/distributors/