



Spinsolve Autosampler




- Easy to mount
- Fully controlled by the Spinsolve software
- 20 slots for 5 mm NMR tubes
- Easy configuration of all protocols for sample analysis
- Compatible with all Spinsolve models
- Affordable and robust



Spinsolve

- Protocols for each sample are added to the queue with a single click
- New samples can be added/removed even during a running queue
- A preset series of experiments can be loaded from a previously saved list



- Easy view of the history of completed and future protocols in the queue
- The sample queue can be modified at any time to change the order in which the samples are measured
- The order of protocols can be easily changed by dragging them up or down

Data is displayed for running and completed experiments

The individual parameters for each experiment are defined when the protocols are loaded and they are displayed during data acquisition

Contact us now for a quote, to request a demo or to measure your samples

Email: sales@magritek.com

Website: www.magritek.com/contact-us

GERMANY +49 241 9278 7270

UNITED STATES +1 855 667 6835

UNITED KINGDOM +44 7468 529 615

NEW ZEALAND +64 4 477 7096

For a complete listing of our global offices and distributors visit: www.magritek.com/about-us/distributors/